

ECE 513, University of Arizona

Heart Rate Monitor System

Project Documentation / User Manual

By
Amber Parker, Seven Gilbert, Samantha Perry

Table of Contents

Table of Contents.....	1
1. Project Overview.....	3
1.1 Project Description.....	3
1.1.1 Intended Users.....	3
1.1.2 System Summary.....	3
1.1.2.1 Frontend (React).....	3
1.1.2.2 Backend (Node / Express / MongoDB).....	3
1.1.2.3 Embedded Device (Photon).....	3
2. System Architecture and Implementation.....	4
2.1 System Architecture Overview.....	4
2.2 Backend Implementation.....	4
2.3 Frontend Implementation.....	5
2.4 Embedded Device Implementation.....	6
3. File and Code Structure Description.....	7
3.1 Backend Files.....	7
3.1.1 Models.....	7
3.1.2 Routes.....	7
3.1.3 Controllers / Services (if applicable).....	8
3.2 Frontend Files.....	8
3.2.1 Public Files.....	8
3.2.2 JavaScript / Components.....	8
3.3 Embedded Device Files.....	9
4. Results and Demonstration.....	9
4.1 Frontend Results.....	10
4.1.1 Login and registration pages.....	10
4.1.2 Patient dashboard.....	13
4.1.3 Heart rate and SpO ₂ visualization.....	15
4.1.4 Device management (add/remove device).....	16
4.1.5 Physician management.....	21
4.2 Embedded Device Results.....	24
4.2.1 Collected data.....	24
4.2.2 Device in operation.....	27
4.2.3 Data correspondence.....	31
4.3 System Workflow Demonstration.....	33
4.3.1 Device boot and configuration.....	33
4.3.2 Measurement capture.....	33
4.3.3 Sending data to the backend.....	36
4.3.4 Frontend consumption.....	38

4.3.5 Physician override.....	38
4.3.6 Offline handling.....	40
5. Lessons Learned.....	41
6. Challenges and Resolutions.....	41
7. Team Contributions.....	42
8. References.....	42
8.1 Framework Documentation.....	42
8.2 Libraries.....	42
8.3 Datasheets.....	43

1. Project Overview

1.1 Project Description

This project is a secure platform designed to monitor heart rate and SpO₂ in a simple and accessible way. Patients use the system to record and view their vital signs, while physicians can monitor their assigned patients and review both daily readings and longer-term weekly trends. Admins and operators are responsible for managing device deployment and overall system setup. The platform also allows measurement frequency to be adjusted as needed, making it flexible for different patient monitoring needs.

The project can be found at this GitHub repository:

<https://github.com/samnperry/ECE513FinalProject>.

1.1.1 Intended Users

Patients can register their device, take measurements, review weekly averages and daily charts, and set measurement frequency to what they want.

Physicians can log into a separate portal to view assigned patients, see 7-day summaries, look at daily readings, and override measurement frequency.

1.1.2 System Summary

1.1.2.1 Frontend (React)

The frontend is a responsive web app built around certain user roles. Patients can create accounts, log in, register their device, and view their heart-rate and SpO₂ data through daily and weekly dashboards, along with basic account management.

Physicians have access to an all-patients view, individual patient summaries, detailed daily readings, and the ability to override how often a patient's device takes measurements when needed.

1.1.2.2 Backend (Node / Express / MongoDB)

The backend provides REST APIs that handle authentication, account updates, device registration, and incoming sensor measurements. It also manages physician-to-patient assignments and holds both patient and physician dashboards. Each device has a unique API key, and physician-defined measurement overrides are stored and sent back to devices through a configuration endpoint (GET /api/device/config/:deviceId).

1.1.2.3 Embedded Device (Photon)

The Photon reads heart-rate and SpO₂ data using the MAX30102 sensor, checks the backend for its current measurement interval, and sends readings to the server using its device-specific API key. If the device loses connectivity, it temporarily stores the data locally and automatically uploads it once the connection is restored.

2. System Architecture and Implementation

2.1 System Architecture Overview

The system is designed as a simple, end-to-end flow that connects a wearable device to a web-based monitoring platform. At a high level, the embedded device (a Particle Photon 2 with a MAX30102 sensor) collects heart rate and SpO₂ data and communicates with a Node/Express backend over HTTP(S). The backend stores and processes this data in MongoDB and makes it available to users through a React-based frontend for both patients and physicians.

The data flow starts on the device, which samples heart rate and SpO₂ values and checks in with the backend to see how often it should be taking measurements. This configuration is fetched from a device-specific endpoint using an API key. Once a valid reading is captured, the device sends the measurement to the backend, again using its API key for authentication. The backend validates the device, stores the data using MongoDB, and exposes it through REST APIs. The frontend then retrieves this information securely and displays it through dashboards and charts tailored to each user role.

The system uses standard web technologies throughout, including HTTPS over Wi-Fi for device communication, JWT-based authentication for users, and a REST-based architecture connecting the backend to the frontend. This approach keeps the system modular, easy to understand, and scalable.

2.2 Backend Implementation

The backend is built using Node.js with the Express framework and MongoDB for data storage, using Mongoose models for users, devices, and measurements. It handles authentication, device management, incoming sensor data, and all patient and physician views. User passwords are hashed with bcrypt, and JWTs are used to secure web-based routes. CORS is configured to allow requests only from the frontend origin.

Authentication routes allow both patients and physicians to sign up and log in, with JWT middleware protecting all authenticated endpoints. Device registration is handled through dedicated device routes, which generate a unique API key for each device. This API key is required for all device communication. Incoming measurements are sent to POST /api/measurements with the device's x-api-key, and each payload includes the device ID, heart rate, SpO₂ value, and timestamp. The backend validates the API key, stores the data in MongoDB, and makes it available through patient and physician dashboard routes.

The backend also manages relationships between patients, physicians, and devices. Patients can assign a physician to their account, and physicians can view summaries for all assigned patients, drill down into daily measurement data, and override how

frequently a device takes measurements. These overrides are stored with the device and returned to the device through the configuration endpoint GET /api/device/config/:deviceId, allowing the embedded system to dynamically update its measurement interval.

Environment variables are used to manage sensitive configuration values, including the database connection string (MONGO_URI), JWT secret (JWT_SECRET), and frontend API base URL (REACT_APP_API_BASE).

For deployment, the React frontend is built locally using npm run build, producing a static build/ directory. This directory is deployed to the AWS EC2 instance under /var/www/heartapp, where it is served by Nginx. The Node.js/Express backend runs on the same EC2 instance on port 5001 and is managed using PM2 to ensure it stays running across logouts and system reboots.

Nginx is configured as a reverse proxy for the application. It serves the static frontend and forwards all /api requests to the backend at 127.0.0.1:5001. HTTPS is enabled using Let's Encrypt via Certbot, which generates and manages SSL certificates for the domain. All traffic is encrypted over port 443, and certificate renewal is handled automatically.

The full deployment was verified using HTTP and HTTPS requests with curl and by reviewing Nginx logs to confirm that authentication and API requests were successfully reaching the backend. This setup provides a complete end-to-end deployment where the frontend is served through Nginx, the backend is reliably managed by PM2, and all communication is secured using HTTPS..

2.3 Frontend Implementation

The frontend is built in React using TypeScript and is structured around role-based access. Patients and physicians see different views and features based on their role, but the overall interface remains consistent and easy to navigate.

Patients can create an account, register or remove devices, assign a physician, and view their data through weekly summary dashboards and detailed daily charts. They also have access to basic account management features, such as updating their password.

Physicians use a separate portal that provides an overview of all assigned patients, including seven-day statistics such as averages, minimums, and maximums. From there, they can view individual patient summaries, drill down into daily measurement details, and adjust the measurement frequency for specific devices when needed.

The frontend communicates with the backend using secure HTTPS requests and includes JWTs in the Authorization header. Data visualization is handled with Chart.js, using line charts for weekly trends and bar charts for daily details, along with tables for

clear presentation of statistics and lists.

- Tech stack/UI: “React + TypeScript with React Router for routing, localStorage for tokens/role, and Pages.css for consistent cards/forms across login/signup/dashboard. Protected routes redirect unauthenticated users.”
- Role routing: “Header/nav hides patient-only links for physicians and vice versa; login redirects physicians to /physician-dashboard and patients to /dashboard.”
- Data flow: “All API calls use fetch with Authorization: Bearer <JWT> and the deployed base https://sfwe513.publicvm.com/api; errors surface as inline banners.”
- Components: “Dashboard cards show weekly stats, measurement frequency controls, and device filters; forms use shared styles; physician tables list patients with action links to summary/daily views; device and physician selectors are populated from live API calls.”
- Charts: “Chart.js line charts for weekly trends, daily detail overlays heart rate vs. SpO2; tables provide min/max/avg and counts for the last 7 days; daily view charts update with date filters.”

2.4 Embedded Device Implementation

The system uses a Particle Photon 2 microcontroller using the Particle Device OS, (Particle.h), with built-in Wi-Fi and two onboard LED's used for visual output indicators. The microcontroller connects to a MAX30102 Heart-Rate / SpO2 sensor via I2C. The sensor is also powered using the voltage output of the Photon 2. The EEPROM hardware storage on the Photon 2 is used for offline data collection.

To collect data, the particle board samples the MAX30102 sensor to read current IR values during the measurement period. When the IR value hits a threshold that determines the presence of a finger, the Photon 2 will start recording data from the sensor with help from the spo2_algorithm.h file included. Data is filled into a rolling buffer that, when populated, is used to calculate the users heart rate and SpO2. These measurements are then recorded and verified by the microcontroller. Once a sufficient number of valid measurements have been taken, the Photon 2 can package and send the data.

The sensor is set to a sampling rate of 25 Hz, or every 40 ms. This data is stored into a rolling buffer of 100 recordings. The data is then used to calculate the heart rate and SpO2 values. Additionally, variables for if those calculations represent valid outputs are set. These need to remain true for a specified number of sequential measures for data to be deemed safe for sending. Time stamps are used for data and event tracking; unix seconds are used to represent our time stamps.

The data is transmitted via WiFi to the Particle Cloud Server. This is done when the custom webhook made in the online Particle Console is triggered using the private particle event, “Photon2_SendEvent”. The microcontroller then uses webhooks to verify whether or not the data was sent successfully. Similarly, the Photon 2 sends a GET request to the webhook using the event “Photon2_Config_Request” in order to obtain the most recent configuration for frequency of measurement times. This information is collected by the front end and saved in the server.

Local processing and filtering is performed on sensor data and recorded measurements. A threshold is set for the incoming IR values to filter out any irrelevant data. Once the sensor begins to make valid measurements, the program waits until the data has stabilized before sending. Measurements made are stored into local memory for up to 24 hours where they can be later uploaded when able.

3. File and Code Structure Description

3.1 Backend Files

3.1.1 Models

Path: server/models

The user.js model stores all user-related information, including email, hashed password, role (patient or physician), registered devices, and any assigned physician. The device.js model represents each registered device and includes the device ID, the associated user, a unique API key, an optional measurement frequency override, and a nickname for easier identification. The measurement.js model stores incoming heart rate and SpO₂ readings along with the device ID and timestamp.

- user.js: Stores users (email, password hash, role: patient/physician, devices[], assignedPhysician).
- device.js: Stores devices (deviceId, user ref, apiKey, measurementFrequencySeconds, nickname).
- measurement.js: Stores heart rate/SpO₂ readings (deviceId, heartRate, spo2, timestamp).

3.1.2 Routes

Path: server/routes

Authentication and account-related functionality is handled in authRoutes.js, which supports signup, login, password updates, and JWT generation. Device-related actions are handled in deviceRoutes.js, including device registration (which returns a unique API key), listing devices, and providing device configuration through /api/device/config/:deviceId. Incoming sensor data is handled by

measurementRoutes.js, which accepts POST requests from devices using an API key and allows recent measurements to be retrieved by device. Physician-specific functionality lives in physicianRoutes.js, covering physician registration and login, patient assignment, all-patient statistics, per-patient summaries, daily views, and device frequency overrides. An optional accountsRoutes.js file is used for additional account-related helpers when needed. Most of the business logic is implemented directly inside the route files. This includes password validation, API key generation during device registration, measurement validation and storage, physician-only data aggregation, saving frequency overrides, and serving configuration data back.

- authRoutes.js: Signup/login, account update (password), JWT issuance.
- deviceRoutes.js: Register device (returns per-device apiKey), list devices, get device config /api/device/config/:deviceId (includes measurementFrequencySeconds).
- measurementRoutes.js: POST measurements with x-api-key, GET recent measurements by device (enforces apiKey/device match).
- physicianRoutes.js: Physician register/login, list physicians, assign physician to patient, all-patients 7-day stats, patient summary/daily, override device frequency.

3.1.3 Controllers / Services (if applicable)

The backend does not use separate controller or service layers. Instead, all business logic is handled directly inside the route files, which keeps the project easier to follow and debug. This includes authentication and password validation, device registration with automatic API key generation, validation, and storage of incoming measurements, and aggregation of patient statistics for physician views. Physician-only functionality, such as overriding device measurement frequency, is also handled here, along with serving updated configuration data back to devices when they request it.

3.2 Frontend Files

3.2.1 Public Files

The public/index.html file serves as the entry point for the React single-page application and includes the landing page introduction shown before the app loads. Static assets such as the favicon, manifest, robots.txt, and team photos are stored directly in the public/ directory and are loaded by the browser as needed.

3.2.2 JavaScript / Components

Utilizing React app is initialized in src/index.tsx, with App.tsx defining the main layout and all route definitions. Authentication components handle patient and physician login and signup. Patient-facing components include the main dashboard, which shows

weekly summaries and daily charts, along with the home page and device registration flow. The account page allows patients to update their password, add or remove devices, and assign a physician. Physician functionality is grouped into its own set of components, including an all-patients dashboard with seven-day statistics, individual patient summaries, detailed daily views, and controls to adjust how often a patient's device takes measurements. Shared components such as the header and navigation bar handle role-based navigation so each user only sees the pages relevant to them. State management is handled mostly with local useState and useEffect hooks. All API communication uses standard fetch calls to the backend defined by REACT_APP_API_BASE, with JWTs passed in the Authorization: Bearer header. Device API keys are never exposed to or used by the frontend.

3.3 Embedded Device Files

The embedded device code is centered around the main firmware file, 513Photon2.ino (with a matching .cpp file generated to contain the same logic). This firmware implements a state machine that controls user prompting, sensor sampling, and data transmission. Device-specific values such as the API key and backend host information are set before building the firmware.

Sensor handling is implemented using the MAX30102 driver and supporting algorithm files, which handle heart rate and SpO₂ calculations. MAX30105.cpp is the SparkFun MAX3010x driver that wraps all I²C register definitions, configuration, FIFO access, LED pulse parameters, etc. It gives the firmware low-level control over the IR/Red photodiode hardware. The file heartRate.cpp implements a Maxim's PBA heart-beat detector: it filters incoming IR samples, tracks AC/DC components, and emits beat events. The spo2_algoryihm.cpp file provides the Maxim SpO₂/heart-rate computation routine (maxim_heart_rate_and_oxygen_saturation) that the firmware calls after filling its sample buffers.

The device uses HTTP networking to request its current configuration from the backend and update its measurement interval dynamically; frequency-webhook.json. Measurements are sent to the backend over HTTP, and the device includes local EEPROM buffering so readings are not lost during temporary network outages. Once connectivity is restored, any queued measurements are automatically sent. 513Photon2.ino handles the implementation of these actions.

4. Results and Demonstration

The system was demonstrated through screenshots showing key functionality for both patient and physician roles. These include login and registration pages, the patient dashboard with weekly summary statistics and daily detail charts, and heart rate and SpO₂ visualizations created using [Chart.js](#). Additional screenshots show device


management from the patient account page, including adding and removing devices, as well as physician management features. The physician portal screenshots highlight the all-patients overview with seven-day averages, minimums, and maximums, individual patient summaries, detailed daily views, and the ability to override device measurement frequency.

4.1 Frontend Results

4.1.1 Login and registration pages

The frontend authentication flow was validated using both the login and registration pages. The login page shows standard email and password inputs along with clear error messaging when invalid credentials are entered. The registration page enforces password requirements and displays a success message once an account is created, confirming that the authentication flow and form validation are working as expected.

Invalid Login:



Valid Login that directs to Dashboard:

Heart Rate Monitor

Dashboard Reference Account Devices

WEEKLY SUMMARY

Dashboard

Pull the last 50 measurements for your device and get a seven-day snapshot of heart rate and SpO2 trends.

Device ID: 0a10aced20219494a064eec Refresh

Register device Assign physician

Stores the device and fetches the latest 50 readings.

TIME RANGE

Filter measurements

Optional start/end times to focus charts on a specific window.

Start: mm/dd/yyyy, --:-- -- End: mm/dd/yyyy, --:-- --

If blank, defaults to last 24 hours for the chart and last 7 days for the summary.

Measurement Frequency

Controls how often the device collects and uploads measurements.

Every 10 minutes

Save frequency

Avg heart rate (7d) — Across last 7 days


Avg SpO2 (7d) — Across last 7 days

Readings collected — From last 7 days

Weekly breakdown

Average heart rate and SpO2 per day. Older than 7 days is hidden.

No data for the last 7 days. Try refreshing.



Sign up without strong password:

Heart Rate Monitor

Dashboard Reference Account Devices Login Sign Up

Sign up

user10@example.com


...

Regular user Physician


Sign up

Password must be at least 8 chars with upper, lower, number, and special character.

Already have an account? [Log in here](#)





Sign up with strong password to dashboard (shows failed to load device information as there is no device registered to the user):



The screenshot shows the "Dashboard" view of the "Heart Rate Monitor" application. At the top, there is a navigation bar with links for Dashboard, Reference, Account, and Devices. The "Dashboard" link is highlighted with a blue background. Below the navigation bar is a "WEEKLY SUMMARY" section with the title "Dashboard". A sub-instruction says "Pull the last 50 measurements for your device and get a seven-day snapshot of heart rate and SpO2 trends." On the left, there is a "Device ID" input field with placeholder "Enter device ID" and a "Refresh" button. Below it are buttons for "Register device" and "Assign physician", with a note that it stores the device and fetches the latest 50 readings. On the right, there is a "TIME RANGE" section titled "Filter measurements" with "Start" and "End" date inputs. A note says "If blank, defaults to last 24 hours for the chart and last 7 days for the summary." Below these sections is a red error message box containing the text "Failed to load device information". Further down are sections for "Measurement Frequency" (set to "Every 30 minutes"), "Avg heart rate (7d)", "Avg SpO2 (7d)", and "Readings collected". The "Readings collected" section shows a value of "From last 7 days". At the bottom, there is a "Weekly breakdown" section with a note about average heart rate and SpO2 per day, stating "No data for the last 7 days. Try refreshing".

Sign up for physician to physician dashboard:



4.1.2 Patient dashboard

After logging in as a patient, the dashboard displays a weekly summary view that includes average heart rate, average SpO₂, and the total number of readings collected. The dashboard also shows when the data was last updated and provides a seven-day

breakdown of measurements for the registered device. This view confirms that data is being correctly pulled from the backend and summarized over the expected time window.

Dashboard:

The screenshot shows the 'Heart Rate Monitor' dashboard interface. At the top, there's a navigation bar with tabs for 'Dashboard', 'Reference', 'Account', and 'Devices'. A 'WEEKLY SUMMARY' section displays the last 50 measurements for a device with ID '0a10aced202194944a064ec'. It includes a 'Refresh' button and a note about fetching the latest 50 readings. Below this is a 'TIME RANGE' section with 'Start' and 'End' date pickers, both set to 'mm/dd/yyyy, --:-- --'. A note says 'If blank, defaults to last 24 hours for the chart and last 7 days for the summary.' A 'Measurement Frequency' section shows 'Every 10 minutes' and a 'Save frequency' button. In the 'Weekly breakdown' section, data is shown for each day from Dec 9 to Dec 15, including average heart rate (e.g., 162 bpm), average SpO2 (e.g., 66%), and the number of readings (e.g., 0 to 5). The 'Daily details (last 24h)' section features a line chart comparing heart rate (blue line) and SpO2 (green line) over a 24-hour period, with both values remaining relatively stable around 160 bpm and 65% respectively.

Backend displaying measurement data:

The screenshot shows the Compass MongoDB interface. On the left, the sidebar displays 'CONNECTIONS (2)' with '402db.w0kxp.mongodb.net' and 'ECE513' selected. Under 'ECE513', there are subfolders: 'admin', 'config', 'hrm', 'devices', 'measurements' (which is highlighted in green), 'users', 'local', and 'test'. The main panel shows the 'measurements' collection with 99 documents. The top navigation bar has tabs for 'measurements' (selected), 'users', and '+'. The sub-navigation bar shows 'ECE513 > hrm > measurements'. Below the tabs are buttons for 'Documents' (99), 'Aggregations', 'Schema', 'Indexes' (1), and 'Validation'. A search bar at the top says 'Type a query: { field: 'value' }'. Below it are buttons for 'Explain', 'Reset', 'Find' (highlighted in green), and 'Options'. The results list five documents, each showing a timestamped measurement record with fields: _id, deviceId, heartRate, spo2, timestamp, and __v.

Document	_id	deviceId	heartRate	spo2	timestamp	__v
1	ObjectId('6940782b19e58f20a83c25eb')	"0a10aced202194944a064eec"	150	76	2025-12-15T21:05:47.712+00:00	0
2	ObjectId('6940783f19e58f20a83c25ee')	"0a10aced202194944a064eec"	150	76	2025-12-15T21:06:07.947+00:00	0
3	ObjectId('6940785319e58f20a83c25f1')	"0a10aced202194944a064eec"	150	76	2025-12-15T21:06:27.560+00:00	0
4	ObjectId('6940786719e58f20a83c25f4')	"0a10aced202194944a064eec"	150	76	2025-12-15T21:06:47.828+00:00	0
5	ObjectId('6940847320fb5cf00150c32b')	"0a10aced202194944a064eec"	214	100	2025-12-15T21:58:11.089+00:00	0

4.1.3 Heart rate and SpO₂ visualization

Daily charts were generated using Chart.js to visualize heart rate and SpO₂ over time for a selected date. The charts clearly show changes throughout the day, with labeled axes for time and measurement values. These visualizations confirm that individual readings are stored correctly and rendered accurately in the frontend.

Charts:

Heart Rate Monitor

Dashboard Reference Account Devices

WEEKLY SUMMARY

Dashboard

Pull the last 50 measurements for your device and get a seven-day snapshot of heart rate and SpO₂ trends.

Device ID

[Refresh](#)

[Register device](#) [Assign physician](#)

Stores the device and fetches the latest 50 readings.

TIME RANGE

Filter measurements

Optional start/end times to focus charts on a specific window.

[Clear](#)

Start: End:

If blank, defaults to last 24 hours for the chart and last 7 days for the summary.

Measurement Frequency

Controls how often the device collects and uploads measurements.

Every 10 minutes

[Save frequency](#)

Avg heart rate (7d)

163 bpm

Across last 7 days

Avg SpO₂ (7d)

81%

Across last 7 days

Readings collected

5

From last 7 days

Weekly breakdown

Average heart rate and SpO₂ per day. Older than 7 days is hidden.

Updated 06:44 PM

DAY	AVG HR	Avg SpO ₂	RANGE (HR)	READINGS
Dec 9	—	—	—	0
WED	—	—	—	0
Dec 10	—	—	—	0
THU	—	—	—	0
Dec 11	—	—	—	0
FRI	—	—	—	0
Dec 12	—	—	—	0
SAT	—	—	—	0
Dec 13	—	—	—	0
SUN	—	—	—	0
MON	163 bpm	81%	150–214 bpm	5
Dec 15	163 bpm	81%	150–214 bpm	5

Daily details (last 24h)

Heart rate and SpO₂ plotted by timestamp to spot trends and spikes.

The graph shows two data series: Heart rate (blue line with circles) and SpO₂ (green line with circles). The x-axis represents time from 02:58 PM to 02:06 PM. The y-axis represents values from 60 to 220. A tooltip at 02:06 PM indicates a heart rate of 150 bpm and an SpO₂ of 76%.

4.1.4 Device management (add/remove device)

The Account page allows patients to register new devices by entering a device ID and nickname, as well as remove devices that are no longer in use. The linked devices list updates immediately after these actions, confirming that device registration and removal are handled correctly by both the frontend and backend.

Add a device (it is added to the device list):

The screenshot shows the 'Account' page with three main sections: Profile, Devices, and Physician. The 'Devices' section contains fields for 'Device ID' (set to 'testdevice') and 'Nickname (optional)' (set to 'test'). A blue 'Add device' button is visible. Below these fields is a note about removal and a list of linked devices: '0a10aced202194944a064eec'. The 'Assign physician' button is present but inactive.

The screenshot shows the 'Account' page with a green header bar indicating 'Device registered'. The 'Devices' section now includes a 'Save changes' button. The 'Device ID' field is empty, and the 'Nickname (optional)' field is also empty. The blue 'Add device' button is still present. The note about removal and the list of linked devices remain the same. The 'Assign physician' button is still present but inactive.

```

_id: ObjectId('693dbfc6db9dc580be50339b')
deviceId : "0a10aced202194944a064eec"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "3cf562803c98eee6f2df540bdc3b45a61e2b93200790ca78dbbfe4c3ce47a38c"
__v : 0
measurementFrequencySeconds : 600

```

```

_id: ObjectId('6940bba3b46a9bb595641a76')
deviceId : "testdevice"
nickname : "test"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "dbb4c25eb9b800c289236430ebf6f233e2dffadb77a4f1d2b928fe699a871fe0"
measurementFrequencySeconds : 1800
__v : 0

```

After refresh of page:

Heart Rate Monitor

Dashboard Reference **Account** Devices

Account

Profile

Update your password. Email changes are not allowed.

Save changes

Devices

Add/remove linked devices and set nicknames.

Add device

Select a device to remove

Note: removal expects the device record ID (from the database). Registering uses the public device ID.

Linked devices:
0a10aced202194944a064eec, testdevice (test)

Physician

Choose a physician to link with your account.

Assign physician

Device removal page:

Heart Rate Monitor

Dashboard Reference **Account** Devices

Account

Profile

Update your password. Email changes are not allowed.

Save changes

Devices

Add/remove linked devices and set nicknames.

Add device

Select a device to remove

Remove device

Note: removal expects the device record ID (from the database).
Registering uses the public device ID.
Linked devices:
0a10aced202194944a064eec,
testdevice (test)

Physician

Choose a physician to link with your account.

Selected physician

Assign physician

Removing a device (It is removed from the list of linked devices):

```
_id: ObjectId('693dbfc6db9dc580be50339b')
deviceId : "0a10aced202194944a064eec"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "3cf562803c98eee6f2df540bdc3b45a61e2b93200790ca78dbbfe4c3ce47a38c"
__v : 0
measurementFrequencySeconds : 600
```

```
_id: ObjectId('6940bacdb46a9bb595641a69')
deviceId : "testdevice"
nickname : "test"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "046f714f65118ff587b4ab34038ed75428b926a0474464a334db4032598ed0de"
measurementFrequencySeconds : 1800
__v : 0
```

Account

Profile

Update your password. Email changes are not allowed.

New password

Save changes

Devices

Add/remove linked devices and set nicknames.

Device ID

Nickname (optional)

Add device

testdevice (test)

Note: removal expects the device record ID (from the database).

Registering uses the public device ID.

Linked devices:
0a10aced202194944a064eec,
testdevice (test)

Physician

Choose a physician to link with your account.

Select physician

Assign physician

Account

Profile

Update your password. Email changes are not allowed.

New password

Save changes

Devices

Add/remove linked devices and set nicknames.

Device ID

Nickname (optional)

Add device

✓ Select a device to remove

0a10aced202194944a064eec

Note: removal expects the device record ID (from the database).

Registering uses the public device ID.

Linked devices:
0a10aced202194944a064eec

Physician

Choose a physician to link with your account.

Select physician

Assign physician

```

_id: ObjectId('693dbfc6db9dc580be50339b')
deviceId : "0a10aced202194944a064eec"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "3cf562803c98eee6f2df540bdc3b45a61e2b93200790ca78dbbfe4c3ce47a38c"
__v : 0
measurementFrequencySeconds : 600


```

4.1.5 Physician management

From the patient account page, a physician can be assigned using a dropdown menu, with a confirmation message displayed after a successful assignment.

Assign a physician:

The screenshot shows the 'Account' page of a web application. At the top, there is a navigation bar with links for 'Dashboard', 'Reference', 'Account' (which is highlighted in blue), and 'Devices'. Below the navigation bar, the page is titled 'Account'. The page is divided into three main sections: 'Profile', 'Devices', and 'Physician'. The 'Profile' section contains a 'New password' input field and a 'Save changes' button. The 'Devices' section contains fields for 'Device ID' and 'Nickname (optional)', and a 'Add device' button. The 'Physician' section contains a 'Select physician' dropdown menu and an 'Assign physician' button. There is also some descriptive text and a note about removing devices.



In the database, the assignedPhysician is the ObjectId to the selected physician:

```
_id: ObjectId('69364cb0ff8c94a67406005e')
email : "testuser@example.com"
passwordHash : "$2b$10$/Xq/awL/ijrYWm5ye0dnre.7iqqj7ssDiEJLgyWZy/nXpyYq0yeem"
role : "user"
▶ devices : Array (1)
  assignedPhysician : ObjectId('69364ccbff8c94a674060061')
  __v : 4
```

```
_id: ObjectId('69364ccbff8c94a674060061')
email : "testphysician@example.com"
passwordHash : "$2b$10$u29ovhrJq4wzTUWv92gR5uDhYICxqBxq8vEviSN9wrP8s1E6vaGYi"
role : "physician"
▶ devices : Array (empty)
  assignedPhysician : null
  __v : 0
```

The physician portal provides an overview of all assigned patients, including seven-day averages, minimums, and maximums for each patient, along with device counts. Physicians can select an individual patient to view a summary of recent readings, access detailed daily charts, and adjust the device measurement frequency when needed. These views confirm that role-based access and physician-specific functionality are working correctly.

Physician dashboard:

The screenshot shows the Physician Dashboard interface. At the top, there is a header with the logo "Heart Rate Monitor", navigation links "Physician Dashboard", "Reference", and "Account", and a share icon. The main content area is divided into two sections: "Physician Dashboard" and "Workflow".

Physician Dashboard:
Backend: <https://sfwe513.publicvm.com>
Manage assigned patients, view summaries, daily details, and override measurement frequency.
[View all patients](#)

Workflow:

- Patients: list assigned patients with device counts.
- Summary: latest readings per device, frequency override controls.
- Daily: pick a date to see time-stamped readings.

All Patients API: <https://sfwe513.publicvm.com/api/physician/patients>

EMAIL	DEVICES	7D AVG	7D MIN	7D MAX	ACTIONS
testuser@example.com	1	158	22	214	Summary Daily

Device Summary and Frequency Adjustment:

The screenshot shows the "Device Summary and Frequency Adjustment" page for the patient "testuser@example.com".

Physician Dashboard:
Backend: <https://sfwe513.publicvm.com>
Manage assigned patients, view summaries, daily details, and override measurement frequency.
[View all patients](#)



Workflow:

- Patients: list assigned patients with device counts.
- Summary: latest readings per device, frequency override controls.
- Daily: pick a date to see time-stamped readings.

testuser@example.com - Device Summary

Patient device: 0a10aced202194944a064eec | Latest: 214 bpm at 2:58:11 PM
Override interval: Frequency updated

Daily view:



4.2 Embedded Device Results

4.2.1 Collected data

Heart rate and SpO₂ readings collected by the device were verified by viewing stored records in the database. Each record includes the device ID, timestamp, heart rate, and

SpO₂ value. The timestamps in the database align with the times reported by the device, confirming that measurements are being transmitted and stored correctly.

Serial monitor of data being read:

```
[MAX30102] IR=420 RED=312 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=34840 RED=12519 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2

[STATE] STATE_ACQUIRE
[MAX30102] IR=36956 RED=32979 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37114 RED=33254 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37186 RED=33527 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37098 RED=33337 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37067 RED=33317 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37062 RED=33335 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37045 RED=33314 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37058 RED=33353 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37040 RED=33330 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37038 RED=33351 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37073 RED=33391 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37032 RED=33363 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37017 RED=33394 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37060 RED=33436 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37044 RED=33445 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37012 RED=33420 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37019 RED=33438 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37035 RED=33491 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37038 RED=33502 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37012 RED=33459 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36974 RED=33445 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36952 RED=33446 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36952 RED=33450 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36961 RED=33490 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36910 RED=33462 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36904 RED=33478 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36904 RED=33486 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36928 RED=33542 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36926 RED=33540 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36902 RED=33509 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36893 RED=33510 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36885 RED=33515 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36865 RED=33546 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3

[STATE] STATE_MEASUREMENT_READY

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
Server ACK received (DB recorded).

[STATE] STATE_FLASH_GREEN
```

```
_id: ObjectId('6940ca0eb46a9bb595641b51')
deviceId : "0a10aced202194944a064eec"
heartRate : 214
spo2 : 100
timestamp : 2025-12-16T02:55:10.513+00:00
__v : 0
```

```
_id: ObjectId('6940ca18b46a9bb595641b54')
deviceId : "0a10aced202194944a064eec"
heartRate : 214
spo2 : 100
timestamp : 2025-12-16T02:55:20.517+00:00
__v : 0
```

```
_id: ObjectId('6940ca77b46a9bb595641b58')
deviceId : "0a10aced202194944a064eec"
heartRate : 42
spo2 : 98
timestamp : 2025-12-16T02:56:55.009+00:00
__v : 0
```

```
_id: ObjectId('6940caecb46a9bb595641b5b')
deviceId : "0a10aced202194944a064eec"
heartRate : 166
spo2 : 39
timestamp : 2025-12-16T02:58:52.180+00:00
__v : 0
```

In addition, the same readings shown in the database appear in the frontend dashboards and charts, demonstrating consistency across the system.

Frontend data:

Heart Rate Monitor

Dashboard Reference Account Devices

WEEKLY SUMMARY

Dashboard

Pull the last 50 measurements for your device and get a seven-day snapshot of heart rate and Sp02 trends.

Device ID

Register device **Assign physician**

Stores the device and fetches the latest 50 readings.

TIME RANGE

Filter measurements

Optional start/end times to focus charts on a specific window.

Clear

Start **End**

If blank, defaults to last 24 hours for the chart and last 7 days for the summary.

Measurement Frequency

Controls how often the device collects and uploads measurements.

Save frequency

Avg heart rate (7d)

159 bpm

Across last 7 days

Avg Sp02 (7d)

84%

Across last 7 days

Readings collected

4

From last 7 days

DAY	AVG HR	AVG SPO2	RANGE (HR)	READINGS
Dec 9	—	—	—	0
WED Dec 10	—	—	—	0
THU Dec 11	—	—	—	0
FRI Dec 12	—	—	—	0
SAT Dec 13	—	—	—	0
SUN Dec 14	—	—	—	0
MON Dec 15	159 bpm	84%	42-214 bpm	4


Daily details (last 24h)

Heart rate and Sp02 plotted by timestamp to spot trends and spikes.


4.2.2 Device in operation

The assembled device consists of a Particle Photon 2 connected to a MAX30102 sensor. LED indicators on the device provide feedback during prompting, measurement, and transmission. Serial monitor output shows the device fetching its configuration, successfully capturing readings, and publishing measurements to the backend.

LED indicating user to place finger:



LED indicating data is being collected:



Serial Monitor of data being collected:

```

[STATE] STATE_PROMPT_USER
[MAX30102] IR=478 RED=333 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=0 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=345 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=467 RED=334 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=460 RED=342 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=461 RED=333 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2
[MAX30102] IR=6616 RED=1144 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=2

[STATE] STATE_ACQUIRE
[MAX30102] IR=37362 RED=33956 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37622 RED=34341 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38163 RED=35068 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38064 RED=34930 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38013 RED=34866 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37997 RED=34810 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37996 RED=34827 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=34823 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37940 RED=34867 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37936 RED=34816 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38035 RED=34875 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38062 RED=34976 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38028 RED=35040 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38044 RED=34978 BP>M=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38042 RED=35045 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38009 RED=34986 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38024 RED=35041 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37970 RED=34972 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38003 RED=34997 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37986 RED=35025 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=35050 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37934 RED=35021 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37921 RED=35015 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37906 RED=34992 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37886 RED=34975 BPM=0 (v=0) Sp02=% (v=0) bufFilled=0 state=3

```

Database with published data:

```
_id: ObjectId('6940ca0eb46a9bb595641b51')
deviceId : "0a10aced202194944a064eec"
heartRate : 214
spo2 : 100
timestamp : 2025-12-16T02:55:10.513+00:00
__v : 0
```

```
_id: ObjectId('6940ca18b46a9bb595641b54')
deviceId : "0a10aced202194944a064eec"
heartRate : 214
spo2 : 100
timestamp : 2025-12-16T02:55:20.517+00:00
__v : 0
```

```
_id: ObjectId('6940ca77b46a9bb595641b58')
deviceId : "0a10aced202194944a064eec"
heartRate : 42
spo2 : 98
timestamp : 2025-12-16T02:56:55.009+00:00
__v : 0
```

```
_id: ObjectId('6940caecb46a9bb595641b5b')
deviceId : "0a10aced202194944a064eec"
heartRate : 166
spo2 : 39
timestamp : 2025-12-16T02:58:52.180+00:00
__v : 0
```

```
_id: ObjectId('6940cc22b46a9bb595641b62')
deviceId : "0a10aced202194944a064eec"
heartRate : 250
spo2 : 92
timestamp : 2025-12-16T03:04:02.774+00:00
__v : 0
```


4.2.3 Data correspondence

Each POST request sent by the device to the /api/measurements endpoint results in a stored record in MongoDB. The device ID and timestamp shown in the serial logs match the corresponding database entry and the data point displayed in the frontend charts. This confirms end-to-end data flow from the device to the user interface.

Database with 100+ collected data:

			100	1 - 100 of 105							
<hr/>											
<pre>_id: ObjectId('693fa6703357bddf4a2bab12') deviceId : "0a10aced202194944a064eec" heartRate : 115 spo2 : 30 timestamp : 2025-12-15T06:10:56.046+00:00 __v : 0</pre>											
<pre>_id: ObjectId('6940782b19e58f20a83c25eb') deviceId : "0a10aced202194944a064eec" heartRate : 150 spo2 : 76 timestamp : 2025-12-15T21:05:47.712+00:00 __v : 0</pre>											
<pre>_id: ObjectId('6940783f19e58f20a83c25ee') deviceId : "0a10aced202194944a064eec" heartRate : 150 spo2 : 76 timestamp : 2025-12-15T21:06:07.947+00:00 __v : 0</pre>											
<pre>_id: ObjectId('6940785319e58f20a83c25f1') deviceId : "0a10aced202194944a064eec" heartRate : 150 spo2 : 76 timestamp : 2025-12-15T21:06:27.560+00:00 __v : 0</pre>											
<pre>_id: ObjectId('6940786719e58f20a83c25f4') deviceId : "0a10aced202194944a064eec" heartRate : 150 spo2 : 76 timestamp : 2025-12-15T21:06:47.828+00:00 __v : 0</pre>											

Chart from frontend:



4.3 System Workflow Demonstration

4.3.1 Device boot and configuration

When powered on, the Photon 2 requests its configuration from the backend using its API key. The device receives the current measurement interval and updates its internal timing accordingly. This update is confirmed through serial log messages indicating the new measurement interval.

4.3.2 Measurement capture

The device prompts the user to place their finger on the sensor and samples data from the MAX30102 at approximately 25 Hz. To ensure accuracy, the firmware requires multiple consecutive stable samples before accepting a measurement. The serial logs show clear state transitions as the device moves from idle, to prompting, to acquiring data, and finally to a ready measurement state.

State transition:

```

Serial monitor opened successfully:
ice ID: 0a10aced202194944a064eec
MAX30102 initialized.

[STATE] STATE_IDLE_WAIT

[STATE] STATE_PROMPT_USER
[MAX30102] IR=478 RED=333 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=0 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=345 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=467 RED=334 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=460 RED=342 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=461 RED=333 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=6616 RED=1144 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2

[STATE] STATE_ACQUIRE

[MAX30102] IR=37362 RED=33956 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37622 RED=34341 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38163 RED=35068 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38064 RED=34930 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38013 RED=34866 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37997 RED=34810 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37996 RED=34827 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=34823 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37940 RED=34867 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37936 RED=34816 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38035 RED=34875 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38062 RED=34976 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38028 RED=35040 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38044 RED=34978 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38042 RED=35045 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38009 RED=34986 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38024 RED=35041 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37970 RED=34972 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38003 RED=34997 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37986 RED=35025 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=35050 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3

```

Multiple stable samples:

```

[STATE] STATE_PROMPT_USER
[MAX30102] IR=478 RED=333 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=0 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=473 RED=345 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=467 RED=334 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=460 RED=342 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=461 RED=333 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2
[MAX30102] IR=6616 RED=1144 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=2

[STATE] STATE_ACQUIRE
[MAX30102] IR=37362 RED=33956 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37622 RED=34341 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38163 RED=35068 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38064 RED=34930 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38013 RED=34866 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37997 RED=34810 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37996 RED=34827 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=34823 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37940 RED=34867 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37936 RED=34816 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38035 RED=34875 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38062 RED=34976 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38028 RED=35040 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38044 RED=34978 BP>M=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38042 RED=35045 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38009 RED=34986 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38024 RED=35041 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37970 RED=34972 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38003 RED=34997 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37986 RED=35025 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37954 RED=35050 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37934 RED=35021 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37921 RED=35015 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37906 RED=34992 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37886 RED=34975 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3

```

Publish measurements:

```

[STATE] STATE_ACQUIRE
[MAX30102] IR=34217 RED=31420 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34355 RED=31612 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34470 RED=31717 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34457 RED=31742 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34596 RED=31894 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34597 RED=31978 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34757 RED=32192 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34691 RED=32078 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34661 RED=32044 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34650 RED=32042 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34667 RED=32053 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34651 RED=32078 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34619 RED=32077 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34634 RED=32068 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34626 RED=32121 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34664 RED=32173 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34631 RED=32125 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34626 RED=32131 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34606 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34596 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34654 RED=32231 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34680 RED=32249 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34616 RED=32206 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34600 RED=32214 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34585 RED=32204 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34593 RED=32211 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34568 RED=32228 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34570 RED=32226 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34611 RED=32240 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34510 RED=32201 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34466 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34488 RED=32165 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34536 RED=32264 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3

[STATE] STATE_MEASUREMENT_READY

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
Server ACK received (DB recorded).

[STATE] STATE_FLASH_GREEN

[STATE] STATE_IDLE_WAIT

```

4.3.3 Sending data to the backend

Once a valid reading is captured, the device sends the measurement to the backend as a JSON payload containing the device ID, heart rate, SpO₂, and timestamp. The backend validates the API key, stores the data in MongoDB, and returns a success response. The device logs confirm receipt of the server acknowledgment.

Payload sent to backend with measurement:

```

○ samanthaperry@Samanthas-MacBook-Pro photon % particle subscribe hook-response/Photon2_SendEvent
t
Subscribing to "hook-response/Photon2_SendEvent" from my devices

{"data":"{\\"deviceID\\":\\"0a10aced202194944a064eec\\",\\"heartRate\\":115,\\"spo2\\":45,\\"_id\\":\\"69
40cf4fb46a9bb595641b66\\",\\"timestamp\\":\\"2025-12-16T03:17:35.300Z\\",\\"_\\"\\":0}","ttl":60,"publ
ished_at": "2025-12-16T03:17:35.322Z", "coreid": "particle-internal", "name": "hook-response/Photon
2_SendEvent/0"}

```

```

bool publishMeasurement(const MeasurementRecord &rec) {
    // Build JSON payload
    String payload = String::format(
        "{\"",
        "\"deviceId\":\"%s\",",
        "\"heartRate\":%d,"
        "\"spo2\":%d,"
        "\"timestamp\":%lu"
        "}",
        deviceId.c_str(),
        (int)rec.heartRate,
        (int)rec.spo2,
        (unsigned long)rec.timestamp
    );

    return Particle.publish(MEAS_EVENT, payload, PRIVATE);
}

```

Device logs with server acknowledgement:

```

[STATE] STATE_ACQUIRE
[MAX30102] IR=34217 RED=31420 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34355 RED=31612 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34470 RED=31717 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34457 RED=31742 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34596 RED=31894 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34597 RED=31978 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34757 RED=32192 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34691 RED=32078 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34661 RED=32044 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34650 RED=32042 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34667 RED=32053 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34651 RED=32078 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34619 RED=32077 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34634 RED=32068 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34626 RED=32121 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34664 RED=32173 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34631 RED=32125 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34626 RED=32131 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34696 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34596 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34654 RED=32231 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34680 RED=32249 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34616 RED=32206 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34600 RED=32214 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34585 RED=32204 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34593 RED=32211 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34568 RED=32228 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34570 RED=32226 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34611 RED=32240 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34510 RED=32201 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34466 RED=32143 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34488 RED=32165 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3
[MAX30102] IR=34536 RED=32264 BPM=187 (v=0) SpO2=81% (v=0) bufFilled=0 state=3

[STATE] STATE_MEASUREMENT_READY

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
Server ACK received (DB recorded).

[STATE] STATE_FLASH_GREEN


[STATE] STATE_IDLE_WAIT

```

4.3.4 Frontend consumption

Both patient and physician views retrieve measurement data through secured REST endpoints using JWT authentication. The dashboards update automatically to reflect new readings, and the time of the posted measurement matches the data point shown in the daily charts.

After data collected from previous section:



4.3.5 Physician override

Physicians can adjust how often a device takes measurements through the frontend interface. When a new frequency is set, the backend stores the updated value and provides it to the device on the next configuration fetch. The device then logs the updated interval and begins using the new measurement cadence.

Original frequency in database:

```
_id: ObjectId('6940c488b46a9bb595641b16')
deviceId : "0a10aced202194944a064eec"
nickname : "test2"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "59613fa10ce3e7d9db64a668216a1aaaf306939fc08e29baf819b6165ac273e26"
measurementFrequencySeconds : 600
__v : 0
```

Physician changes frequency:

The screenshot shows the 'Heart Rate Monitor' application interface. At the top, there is a navigation bar with links for 'Physician Dashboard', 'Reference', 'Account', and a share icon. The main area is divided into two sections: 'Physician Dashboard' on the left and 'Workflow' on the right. The 'Physician Dashboard' section contains a 'Backend' link to 'https://sfwe513.publicvm.com' and a note about managing assigned patients. It includes a button labeled 'View all patients'. The 'Workflow' section lists three items: 'Patients: list assigned patients with device counts.', 'Summary: latest readings per device, frequency override controls.', and 'Daily: pick a date to see time-stamped readings.' Below these sections, a user summary for 'testuser@example.com' is shown, indicating a patient device with ID '0a10aced202194944a064eec' and a latest reading of 115 bpm at 8:17:35 PM. An 'Override interval' dropdown set to '30 min' has an 'Apply' button next to it, with a green status message 'Frequency updated'.

Frequency changes in the backend:

```
_id: ObjectId('6940c488b46a9bb595641b16')
deviceId : "0a10aced202194944a064eec"
nickname : "test2"
user : ObjectId('69364cb0ff8c94a67406005e')
apiKey : "59613fa10ce3e7d9db64a668216a1aaf306939fc08e29baf819b6165ac273e26"
measurementFrequencySeconds : 1800
__v : 0
```

Device recognizing frequency update and sending GET request to retrieve new frequency:

```

[MAX30102] IR=37864 RED=34782 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38042 RED=34812 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=39223 RED=36102 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=39627 RED=36643 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38861 RED=35574 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=39171 RED=36273 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=39052 RED=35577 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38281 RED=34378 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38045 RED=34080 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37706 RED=33765 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3

[STATE] STATE_MEASUREMENT_READY

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
Server ACK received (DB recorded).

[STATE] STATE_FLASH_GREEN

[STATE] STATE_IDLE_WAIT
Config response received: {"measurementFrequencySeconds":1800}
Updated measurement interval: 1800 sec

```

4.3.6 Offline handling

If the device temporarily loses connectivity, measurements are stored locally instead of being discarded. Once the connection is restored, the queued readings are sent to the backend and appear in both the database and the frontend dashboards. Serial logs confirm when data is stored locally and when it is successfully flushed.

Data collection when WiFi is turned off, then turned back on:

```

[MAX30102] IR=36929 RED=33519 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36954 RED=33619 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36931 RED=33648 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=36967 RED=33638 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=39529 RED=33117 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=40257 RED=37503 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=40074 RED=37204 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38751 RED=35903 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37219 RED=34508 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37257 RED=34233 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=37235 RED=34472 BPM=0 (v=0) SpO2=% (v=0) bufFilled=0 state=3
[MAX30102] IR=38967 RED=35485 BPM=150 (v=1) SpO2=21% (v=1) bufFilled=1 state=3

[STATE] STATE_MEASUREMENT_READY

[STATE] STATE_STORE_OFFLINE
Stored measurement offline in EEPROM queue.

[STATE] STATE_FLASH_YELLOW

[STATE] STATE_IDLE_WAIT

[STATE] STATE_FLUSH_BACKLOG

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
HOOK ERROR event=hook-error/Photon2_SendEvent/0 data=queryAaaa ETIMEOUT sfwe513.publicvm.com
Server returned hook-error.

[STATE] STATE_IDLE_WAIT

[STATE] STATE_FLUSH_BACKLOG

[STATE] STATE_SEND_PENDING
Published measurement event (waiting for server ACK)...

[STATE] STATE_WAIT_SERVER_ACK
Server ACK received (DB recorded).

[STATE] STATE_FLASH_GREEN

[STATE] STATE_IDLE_WAIT

```

5. Lessons Learned

- We learned that a small mismatch in REACT_APP_API_BASE (HTTP vs HTTPS, localhost vs domain) leads to mixed-content blocks and failed API calls; rebuilding with the correct base and redeploying the static build fixed the issue.
- We learned how important it is to set up gitignore files super early on to avoid pushing large amounts of unnecessary data especially when it comes to photon files. We made the mistake of pushing some of these files and had to revert back.
- We learned that having a domain name was essential. It allowed us to obtain valid TLS certificates for the web app and gave the Photon a stable host to call without IP changes. Using the domain mapped to the Elastic IP meant we could terminate HTTPS for browsers while keeping port 80 open for the device, avoiding hardcoded IPs and certificate mismatches.
- We learned that serving a React SPA behind Nginx requires proper proxy and try_files setup; incomplete server blocks or Certbot edits led to TLS/proxy errors until we replaced them with full 80/443 configs.
- We learned that PM2 + systemd for the backend and Nginx for static assets is a simple, reliable deployment model on EC2. It keeps services running after logout/reboot.

6. Challenges and Resolutions

- Challenge 1: Mixed-content errors when the HTTPS page called HTTP APIs.
Resolution: Hardcoded the API base to <https://sfwe513.publicvm.com> and rebuilt the frontend.
- Challenge 2: Certbot left partial Nginx blocks causing ERR_SSL_PROTOCOL_ERROR.
Resolution: Replaced with complete 80/443 server blocks including root/proxy and cert paths, then reloaded Nginx.
- Challenge 3: Backend stopping after SSH logout.
Resolution: Ran PM2 with pm2 startup and pm2 save to auto-start and persist the Node process.
- Challenge 4: Photon needed a stable host and per-device API key.
Resolution: Set API_HOST to the domain after Nginx was set up, kept port 80 open, added /api/device/config/:deviceId for frequency, and enforced apiKey/device checks in ingestion routes.
- Challenge 5: Photon attempted to request configuration before establishing a

cloud connection, resulting in missed webhook triggers. Resolution: Added a cloud-connection check and delayed configuration requests until Particle.connected() returned true.

7. Team Contributions

Team Member	Frontend (%)	Backend (%)	Embedded Device (%)	Documentation (%)	Demos/Testing (%)
Amber	45	10	5	25	15
Seven	5	10	55	15	15
Samantha	10	45	10	15	20

8. References

8.1 Framework Documentation

- [React Framework](#)

8.2 Libraries

- [heartRate.cpp](#)
- [MAX30105.cpp](#)
- [spo2_algorithm.cpp](#)
- React – Frontend JavaScript library for building user interfaces
- React DOM – DOM rendering for React applications
- React Router DOM – Client-side routing for React
- React Scripts – Build and development tooling for React applications
- TypeScript – Typed superset of JavaScript used across the project
- Express.js – Backend web framework for building REST APIs
- Mongoose – Object Data Modeling (ODM) library for MongoDB
- MongoDB – NoSQL database used for storing user and device data
- CORS – Middleware for enabling cross-origin requests
- bcrypt – Library for hashing and securing user passwords
- jwt-simple – Lightweight JSON Web Token (JWT) authentication library
- Material UI (MUI) – React component library for UI design
- MUI Icons – Icon set for Material UI components
- Emotion – CSS-in-JS styling library used with Material UI

- Chart.js – Library for rendering charts and visualizing heart-rate data
- Testing Library (React, DOM, User Event) – Tools for testing React components and user interactions
- Web Vitals – Library for measuring frontend performance metrics

8.3 Datasheets

- [Particle Photon 2 Datasheet](#)
- [MAX30102 Sensor Datasheet](#)